
GAS User Guide


Ahmed Karam


ahmedkaram.edu@gmail.com

Tips


- If the first attempt was unsuccessful (i.e. the first step in the blast took a long time),
 try again after some time.
- GAS deals with NCBI, so in the event of any interruption or delay in the process, this
 is caused by NCBI and not by GAS.
- Use GAS with blast instead of NCBI blast online for large genomes and for assemblies defined by accessions. So, if the size of the genome you are working on is 12 Mb and you have 100 assemblies, the NCBI blast online can do the process without errors, and then use GAS to purify the results.
- Make sure that the "size" folder is in the same folder as GAS.py because it depends on it. That is, when you download the package, do not try to change their location.








The first option to enter the data is to enter the scientific name and it must be correct. Copy correct scientific name and paste (or write it correctly).




The second option for entering data is assembly accession or accessions. Search by scientific name in the assembly database, choose the assemblies, download the ID table, and browse the table through GAS. (you can merge a set of ID tables)



Enter your email and this is what NCBI requires to contact you if anything happens.




Browse a FASTA file that contains sequences




It filters all XML files that are similar to the name of the entered FASTA file that are in the same folder without blastn action.



Choose megablast or blastn





Follow the steps of the GAS through the CMD or the terminal